

University of Bristol

DNA SAMPLE for Angela Hibbert's research programme

OWNER

Name:
Address:
Phone No:
Email:

CAT INFORMATION

Name:
Registration No:
Microchip:
Colour:
Date of Birth:
Sex:

To take a swab sample you will need cotton buds with cotton wool both ends

- 1 Make sure the cat has not eaten or drunk for 30 minutes before taking the mouth swab.
- 2 Rub and twist the cotton bud along the inside of your cats' cheek for about 10 seconds.
- 3 Air dry the swab for about 30 seconds and place in a plastic zip lock bag.
- 4 Repeat the above process using a second cotton bud and place in the zip lock plastic bag.
- 5 Write the cats name clearly on the bag or label the front of each bag
- 6 Complete the form.
- 7 Tape the bag to the form along with a copy of your cat's pedigree.
- 8 If you are going to swab more than one cat please complete the process for each cat separately.
- 9 Send all the paperwork to: